



TANNIN KIT

OCTA-SLIDE 2, 1-10 ppm

CODE 7831-01


QUANTITY	CONTENTS	CODE
15 mL	*Tannin Reagent #1	*7833-E
60 mL	*Tannin Reagent #2	*7834-H
2	Test Tubes, 2.5-10 mL, plastic, w/caps	0106
1	Pipet, 1 mL, plastic	0354
1	Tannic Acid Octa-Slide 2 Bar, 1-10 ppm	7832-01
1	Octa-Slide 2 Viewer	1101

*WARNING: Reagents marked with an * are considered to be potential health hazards. To view or print a Material Safety Data Sheet (MSDS) for these reagents go to www.lamotte.com. To obtain a printed copy, contact LaMotte by e-mail, phone or fax.

To order individual reagents or test kit components, use the specified code number.

The color standards in the LaMotte Tannin Test Kit are calibrated in terms of parts per million tannic acid. Lignin produces about one-third the amount of color that is given by tannin or tannic acid. The test is not specific for tannin or lignin and unless tannin or lignin is definitely known to be present the results may be reported in more general terms such as “tannin-like” or “lignin-like” or simply “hydroxylated aromatic” compounds.








USE OF THE OCTA-SLIDE 2 VIEWER



The Octa-Slide 2 Viewer should be held so non-direct light enters through the back of the Viewer. Insert the reacted sample into the top of the Viewer. Slide the Octa-Slide 2 Bar into the Viewer and match the color of the reaction to the color standards.

WARNING! This set contains chemicals that may be harmful if misused. Read cautions on individual containers carefully. Not to be used by children except under adult supervision

PROCEDURE

<p>1.</p>  <p>Fill a test tube (0106) to 5 mL line with sample water.</p>	<p>2.</p>  <p>Add 3 drops of *Tannin Reagent #1 (7833). Cap and mix.</p>
<p>3.</p>  <p>Use the 1 mL pipet (0354) to add 1 mL of *Tannin Reagent #2 (7834). Cap and mix.</p>	<p>4.</p>  <p>Allow to stand for 30 minutes.</p>
<p>5.</p>  <p>Insert the Tannic Acid Octa-Slide 2 Bar (7832-01) into the Octa-Slide 2 Viewer (1101).</p>	<p>6.</p>  <p>Insert test tube into the Octa-Slide 2 Viewer.</p>
<p>7.</p>  <p>Match sample color to a color standard. Record as ppm Tannic Acid.</p>	

LaMOTTE COMPANY

Helping People Solve Analytical Challenges®

PO Box 329 • Chestertown • Maryland • 21620 • USA
800-344-3100 • 410-778-3100 (Outside U.S.A.) • Fax 410-778-6394

Visit us on the web at www.lamotte.com