



A New Spin on Precision Water Analysis





See Video at
www.lamotte.com/spinvideo



WaterLink® Spin Lab

Code 3576

Let the innovative WaterLink® Spin photometer do all your pool and spa water testing for you. Just fill one unique Spin reagent disk with water and vital tests are done automatically [see chart].

In just 60 seconds all the test results are transferred into the DataMate software (sold separately) and analyzed to display results and recommend precise treatment instructions. Remove the disk and you are ready for the next customer!

This ground-breaking analysis system is so simple anyone can use it! Each sealed reagent disk contains the precise amount of reagent needed to run a complete series of tests. No measuring water, no vials to fill, no prep time!

- So Simple Anyone Can Use It
- Pre-Measured Reagents
- Just Fill Disk With Pool or Spa Water and Spin
- All Test Results in Just 60 Seconds
- Precise Wet Chemistry



Chlorine/Bromine Series Disk

Code 4328-H [50/pk]

Test Factor	Range
Free Chlorine	0-15.0 ppm
Total Chlorine	0-15.0 ppm
Bromine	0-33.0 ppm
pH	6.3-8.6 pH
Calcium Hardness	0-1200 ppm
Total Alkalinity	0-250 ppm
Cyanuric Acid	5-150 ppm
Copper	0-3.0 ppm
Iron	0-3.0 ppm

Includes 100 Chlorine disks with lab

Chlorine/Bromine plus Borate Series Disk

Code 4330-H [50/pk]

Test Factor	Range
Free Chlorine	0-15.0 ppm
Total Chlorine	0-15.0 ppm
Bromine	0-33.0 ppm
pH	6.3-8.6 pH
Calcium Hardness	0-1200 ppm
Total Alkalinity	0-250 ppm
Cyanuric Acid	5-150 ppm
Copper	0-3.0 ppm
Iron	0-3.0 ppm
Borate	0-60 ppm

Disks sold separately

Biguanide plus Borate Series Disk

Code 4331-H [50/pk]

Test Factor	Range
Biguanide	0-70 ppm
Biguanide Shock	0-250 ppm
pH	6.3-8.6 pH
Calcium Hardness	0-1200 ppm
Total Alkalinity	0-250 ppm
Copper	0-3.0 ppm
Iron	0-3.0 ppm
Borate	0-60 ppm

Disks sold separately

WaterLink® DataMate Web

WaterLink® DataMate Web is a cloud-based pool and spa water analysis system that can be used in the store or on-site.

Customize

The program generates custom pool and spa water treatment recommendations that are listed in the proper treatment sequences with precise dosages and clear instructions. These can be viewed on your PC, smartphone or tablet device. An unlimited number of pool and spa tests are saved with the customer record by date and time for easy access and reference.

Data


More control and flexibility than ever, DataMate Web offers pool and spa professionals an incredible level of customization. Add unlimited products to the sanitizers, shocks, balancers, algaecides and specialty chemicals sections—even customize the ideal test factor levels. An equipment log for each customer stores information pertaining to manufacturer, brand, serial number, installation date and warranty expiration of equipment items like heaters, pumps, motors and many more.

Link

The DataMate Web program is compatible with our WaterLink® Spin, Bluetooth® Mobile WaterLink® Spin and WaterLink® 3 Express Lab.

Analyze

The test entry section of the DataMate Web includes fields for Free Chlorine, Total Chlorine, Bromine, Biguanide, Biguanide Shock, pH, Alkalinity, Calcium Hardness, Cyanuric Acid, Copper, Iron, Nitrate, Borate, Phosphate, Manganese, Salt and TDS. There are check boxes for certain water conditions such as water color, stains, cloudy, foamy, or products needed to remedy problems like enzymes, cover cleaner, filter cleaner, surface cleaner, pool opening and closing procedures, weekly maintenance and salt maintenance procedures.



Software Requirements:

For use with WaterLink® DataMate Web from a PC requires the use of Google Chrome™ and Microsoft .NET® Framework 4.0.

Hardware Requirements for use with a WaterLink® Spin Lab:

Recommended Minimum: 3 GHz Intel® Pentium® IV, Core 2 processor or higher with 1 GB RAM or more. Minimum disk space: x86 [32 bit] – 850 MB x64 [64 bit] – 2 GB. High Speed USB 2.0 Port Available.


Internet Explorer® and Microsoft .NET® are registered trademarks of the Microsoft® group of companies.

Google Chrome™ is a registered trademark of Google Inc.

Bluetooth® is a registered trademark of Bluetooth® SIG, Inc.



Go to www.waterlinkconnect.com to Learn More



DataMate 10 Software

Code 1768 [Compact Disk]

- Results Sent Directly From **WaterLink® Spin**
- Easily Customized by Retailer
- Equipment Log

The DataMate 10 software is a user-friendly Windows® PC-based water analysis computer program that generates custom pool and spa water treatment recommendations for your customers. Recommendations are in proper sequence with precise dosages and clear instructions. The program allows more control and flexibility for the pool and spa professional to make the necessary changes to fit their needs.

Software Requirements:

Supported Operating Systems: Microsoft® Windows® 7; Windows Server® 2003 Service Pack 2; Windows Vista® Service Pack 1; Windows® XP Service Pack 3. Prerequisites: Windows® Installer 4.5 or later; Internet Explorer® 5.01 or later.

Hardware Requirements:

Recommended Minimum: 3 GHz Intel® Pentium® IV, Core 2 processor or higher with 1 GB RAM or more. Minimum disk space: x86 [32 bit] – 850 MB x64 [64 bit] – 2 GB. High Speed USB 2.0 Port Available.

Windows®, Windows® XP, Windows Vista®, Windows® 7, Internet Explorer®, Windows Server® and SQL Server®, Microsoft .NET®, Microsoft Excel®, and Microsoft Access® are registered trademarks of the Microsoft® group of companies.

Pentium® is a registered trademark of the Intel® corporation in the U.S. and other countries.

Customize

Insert three pool chemical lines and two spa chemical lines to reflect your brands. Recommendations for sanitizers, shocks, balancers, algaecides and specialty chemicals can be edited along with ideal test factor levels. A new feature is the ability for the retailer to add their logo to the test report.

Data

Store data on your customers like the brand of chlorine [salt] generators, alternative sanitizers, auto cleaners, water features, filters, pool surface, spa model, heater, motor and many other fields. An equipment log for each customer stores information pertaining to manufacturer, brand, serial number, installation date and warranty expiration. Extra fields allow the retailer to add comments on the customer's report or private comments in their file.

Link

The DataMate 10 software is compatible with the WaterLink® Spin, WaterLink®3 Express and WaterLink® Express. Test results from the meter can be sent directly into the software through a USB connector cable or results can be manually typed into the software for tests such as Salt, Phosphates or TDS levels.

Analyze

The test entry section of the DataMate 10 includes fields for Free Chlorine, Total Chlorine, Bromine, Biguanide, Biguanide Shock, pH, Alkalinity, Calcium Hardness, Cyanuric Acid, Copper, Iron, Nitrate, Borate, Phosphate, Manganese, Salt and TDS. There are check boxes for certain water conditions such as water color, stains, cloudy, foamy, or products needed to remedy problems like enzymes, cover cleaner, filter cleaner, surface cleaner, pool opening procedures, pool closing procedures, weekly maintenance and salt maintenance procedures.

Storage

Unlimited pool and spa tests can be stored in customer's files for easy access and reference. In addition, the database can sort information, generate reports, graph specific test factor trends and print labels.

All printed recommendations feature the dealer's store information—name, address, phone & fax numbers, website, etc., as well as a section to advertise monthly sales or other pertinent information for specific customers.



Water Testing Leader Since 1919!

LaMotte Company · PO Box 329 · Chestertown · Maryland · 21620 · USA
t: 800-344-3100 · 410-778-3100 · f: 410-778-6394 · www.lamotte.com/spin